

Special Offer for Academic Institutions

The Hitachi LaChrom Elite® HPLC System with Academic Discount AND Free Software

For a Limited Time Only, with a purchase of a new Hitachi LaChrom Elite® HPLC system, you can receive Hitachi's new **D-2000 Elite Chromatography Software AT NO EXTRA COST** in addition to an **academic discount!***

LaChrom Elite® High-Performance Liquid Chromatography System

Precision, Flexibility and Reliability for the Most Demanding HPLC Applications


Designed to meet the challenges of today's laboratory environment, fast results and demanding workload, LaChrom Elite® couples automated tasks with innovative design.

- ✓ **Flexibility**
Accommodates a wide variety of conventional HPLC applications.
- ✓ **Innovation**
Pump-autosampler synchronization, built-in degasser, programmable wavelength switching and more.
- ✓ **Performance**
Robustness, reliability, and superior performance in your laboratory for many years to come.



Introducing Hitachi D-2000 Elite Chromatography Data System

High-Value Software for Hitachi's LaChrom Elite® HPLC




Example of results transferred to Excel®

Automated Exporting of Data to Excel®

A highly beneficial feature of the D-2000 Elite software is the ability to automatically export report information into Excel using the built-in macros.

2-Peak Deconvolution Functionality

The deconvolution feature mathematically "separates" two overlapping peaks for effective qualitative and quantitative analysis.



- ✓ **Simple User Interface**
Based on the L-7000 HPLC System Manager software, the D-2000 Elite has a simple user interface that is easy to use.
- ✓ **LaChrom Elite® Control**
The D-2000 Elite is a cost-effective software alternative for controlling up to two LaChrom Elite® HPLC systems on a single PC.
- ✓ **3-D DAD Processing Included**
Full control of the L-2455 Diode Array Detector with 3-D data processing is standard. No software add-ons required.
- ✓ **Flexible Report Layouts**
Reports are configured to meet your needs using the flexible reporting functionalities.

Hitachi—The Value Leader in Liquid Chromatography

Hitachi High Technologies America, Inc.
Life Sciences Division
5100 Franklin Drive, Pleasanton, CA 94588

Toll Free: (800) 548-9001
E-mail: sales-LS@hitachi-hta.com
Web site: www.hitachi-hta.com/LSHome
Online store: store.hitachi-hta.com

*Offer valid in the U.S. and cannot be combined with any other offer or discount. Contact Hitachi for details. Requires the purchase of a minimum of three LaChrom Elite HPLC modules (pump, detector plus another). Services excluded from academic discount. The D-2000 software requires Windows XP SP2. Control of L-2455 DAD is limited to one per computer. Hitachi is a registered trademark of Hitachi, Ltd. LaChrom Elite is a registered trademark of Hitachi High Technologies America, Inc. Excel is a registered trademark of Microsoft Corporation. Copyright © 2009 Hitachi High Technologies America, Inc. All rights reserved.